



"Your R&D Partner"

## DATA SHEET

<b>Product Name:</b>	15N Beta-Amyloid (1-38)
<b>Catalog #:</b>	A-1147-1
<b>Source:</b>	Recombinant. A DNA sequence encoding the human beta-amyloid (1-38) sequence was expressed in E. Coli with N15 as the nitrogen source.
<b>Sequence:</b>	<b>D A E F R H D S G Y E V H H Q K L V F F A</b> Asp-Ala-Glu-Phe-Arg-His-Asp-Ser-Gly-Tyr-Glu-Val-His-His-Gln-Lys-Leu-Val-Phe-Phe-Ala- <b>E D V G S N K G A I I G L M V G G</b> Glu-Asp-Val-Gly-Ser-Asn-Lys-Gly-Ala-Ile-Ile-Gly-Leu-Met-Val-Gly-Gly
<b>Molecular Weight:</b>	4181.7 Da by Mass Spec Analysis (theoretical 4181.64)
<b>Peptide Purity:</b>	>97%
<b>Peptide Content:</b>	0.1 mg/vial; (actual peptide content = 0.1 mg)
<b>Counter Ion:</b>	TFA
<b>Supplied As:</b>	White lyophilized powder
<b>Resuspension:</b>	Resuspend in 1 % NH <sub>4</sub> OH, at a concentration of 1 mg/ml. Sonicate for 30 seconds to 1 minute after it has gone into solution. <u>To bring it into your buffer:</u> After resuspension, add 5x or 10x buffer stock and water, to bring to 1x buffer.
<b>Storage:</b>	-20°C
<b>Description:</b>	Beta-amyloid peptide (Aβ), the major constituent of amyloid plaques in the brains of Alzheimer's patients, is thought to be the cause of Alzheimer's Disease (AD). AD is the most common neurodegenerative disease and afflicts about 10% of the population over 60 <sup>4</sup>
<b>References:</b>	1) Yankner, BA, et. al., (1990) Science, <b>250</b> : 279-282 2) Selkoe, D.J., (2001) Physiol. Rev, <b>81</b> : 741-766 3) Stine, W.B. et. al., (2003) J. Biol. Chem, <b>278</b> : 11612-11622 4) Frank, R.A., et. al., (2003) Neurobiology of Aging, <b>24</b> : 521-536

**For research purposes only. Not for human use.**

[www.rpeptide.com](http://www.rpeptide.com)

4241 Mars Hill Road • Bogart • Georgia • 30622 • USA  
Tel 678-753-0747 • Fax 678-753-0746 • E-mail [info@rpeptide.com](mailto:info@rpeptide.com)