



"Your R&D Partner"

DATA SHEET

Product Name: Rabbit Anti-Tau polyclonal antibody
Catalog #: T-1308
Volume/quantity: 0.1mg
Product form: Semipurified antiserum - liquid
Preservatives: 0.1% Thimerosal
Immunogen: recombinant human six isoforms of tau protein
Immunization: Multiple s.c. injections of rabbit
Specificity: Recognizes recombinant and native tau protein
Species reactivity: Human, pig, rat. Others not tested
Mol. Wt. of Antigen: 50 kDa
Positive control: Recombinant human tau protein, brain tissue

Applications and recommended dilutions:

Technique	Yes	Recommended dilutions
ELISA	+	1:40000
Western Blotting	+	1:40000
Immunohistology	+	1:1000

Where this antibody has not been tested for use in a particular techniques this does not necessarily exclude its use in such procedures.

Storage conditions: -20°C (should be stored undiluted). Avoid repeated freezing and thawing.

For research use only. Not for use in humans.

www.rpeptide.com

4241 Mars Hill Road • Bogart • Georgia • 30622 • USA
Tel 678-753-0747 • Fax 678-753-0746 • E-mail info@rpeptide.com