

DATA SHEET

Product Name: Biotin-LC Beta-Amyloid (1-42)

Catalog #: A-1118

Source: Synthetic peptide

Molecular Mass: 4,853.6 Da

Protein Purity: >95%

Counter Ion: TFA

Supplied As: White lyophilized powder

Resuspension: Resuspend in 1% NH₄OH at conc. of .1-1 mg/ml. Recommended to briefly centrifuge to ensure full resuspension of product.

Storage: -20°C

Description:

Biotin-LC (N-terminus) labeled beta-amyloid peptide (Abeta), (*Linker Chain) Labeled (*6 carbon linker - amino hexanoic acid). Beta-amyloid (A-beta) has been long reported as the major constituent of amyloid plaques in the brains of Alzheimer's patients and is believed by many to be the cause of Alzheimer's Disease (AD). AD is the most common neurodegenerative disease and afflicts more than 10% of the population over 65. Recombinantly expressed and sourced from E. coli, rPeptide's high quality beta-amyloid products offer batch-to-batch consistency and ultrapure starting material for your research needs. The TFA (trifluoro-acetic acid) counter-ion is a peptide which was lyophilized in a final solution containing TFA and is a popular choice for researchers wishing to build on previous studies using this same product.

References:

1. Yankner, B., et al., (1990) Science, 250: 279-282
2. Stine, W., et al., (2003) J. Biol. Chem, 278: 11612-11622
3. Frank, R., et al., (2003) Neurobiology of Aging, 24: 521-536
4. Gericke, C., et al., (2020) Front. Immunol., 11: 1029

For research use only. Not for use in humans.



For research use only. Not for use in humans.

www.rpeptide.com

1050 Barber Creek • Bldg: 300 • Suite: 103 • Watkinsville • Georgia • 30677 • USA
Tel 678-753-0747 • Fax 678-753-0746 • E-mail sales@rpeptide.com