

DATA SHEET

Product Name: Alpha Synuclein 10D2 Antibody

Catalog #: SA-3003

Purity: > 95% by SDS-PAGE

Host/Isotype: Mouse / IgG_{1k}

Supplied As: PBS pH 7.4 without additives

Dilution: 2 ng/ml in ELISA, 1-2 µg/ml Western Blot, 1-5 µg/ml in IHC

Storage: -20 °C (repeated thawing and freezing should be avoided)

Description: Alpha Synuclein 10D2 is a monoclonal antibody that is specific for human alpha synuclein. 10D2 has been shown to work within ELISA, Western Blot, as well as Immunohistochemistry experiments.

For research use only. Not for use in humans.