

DATA SHEET

Product Name: Alpha Synuclein pAb 63 Antibody

Catalog #: SA-3004

Product: Polyclonal antibody

Host: Rabbit

Immunogen: Recombinant human Alpha Synuclein

Purity: Undiluted serum

Storage: -20 °C (repeated thawing and freezing should be avoided)

Application: ELISA, Western Blot

Dilution: 1:2.000-10.000 in ELISA, 1:1.000 in Western Blot

Reactivity: Human Alpha Synuclein

For research use only. Not for use in humans.