

DATA SHEET

Product Name: Tau 10F8 Antibody

Catalog #: TA-3007

Purity: > 95% by SDS-PAGE

Host/Isotype: Mouse / IgG2b_k

Supplied As: PBS pH 7.4 without additives

Dilution: 1 µg / ml Western Blot

Storage: at -20 °C (repeated thawing and freezing should be avoided)

Description: Tau 10F8 is a monoclonal antibody that reacts with Tau that is phosphorylated at amino acid 202 (T). Tau 10F8 has been shown to work within ELISA, Western Blot, as well as Immunohistochemistry experiments.

For research use only. Not for use in humans.