

DATA SHEET

Product Name: Tau 2B11 Antibody

Catalog #: TA-3008

Purity: > 95% by SDS-PAGE

Host/Isotype: Mouse / IgG2b_k

Reactivity: Human TAU phosphorylated at amino acid position 231 (T)

Supplied As: PBS pH 7.4 without additives

Dilution: 1 µg / ml Western Blot

Storage: at -20 °C (repeated thawing and freezing should be avoided)

Description: Tau 2B11 is a monoclonal antibody that reacts with Tau that is phosphorylated at amino acid 231 (T). Tau 2B11 has been shown to work within ELISA, Western Blot, as well as Immunohistochemistry experiments.

For research use only. Not for use in humans.