

**DATA SHEET**

**Product Name:** Tau 5G7 Antibody

**Catalog #:** TA-3009

**Purity:** > 95% by SDS-PAGE

**Host/Isotype:** Mouse / IgG2b<sub>k</sub>

**Supplied As:** PBS pH 7.4 without additives

**Dilution:** 1 µg / ml Western Blot

**Storage:** at -20 °C (repeated thawing and freezing should be avoided)

**Description:** Tau 5G7 is a monoclonal antibody that reacts with Tau that is phosphorylated at amino acid 231 (T). Tau 5G7 has been shown to work within ELISA, Western Blot, as well as Immunohistochemistry experiments.

***For research use only. Not for use in humans.***