

DATA SHEET

Product Name: Beta-Amyloid mab 6D11 Antibody

Catalog #: TA-3021

Purity: > 95% by SDS-PAGE

Host/Isotype: Mouse / IgG1_k

Supplied As: PBS pH 7.4 without additives

Dilution: 2 ng/ml in ELISA, 1 µg/ml in Western Blot, 1-5 µg/ml in IHC

Storage: At -20 °C (repeated thawing and freezing should be avoided)

Immunogen: Beta-amyloid (1-42)

Specificity: Beta-amyloid pathological related deposits, reacted with peptide fragment (1-11) of beta-amyloid

Reactivity: Beta-amyloid pathological related deposits

For research use only. Not for use in humans.